SAFE USE OF CRUTCHES BEFORE USE

Please ensure the crutches are adjusted to the correct height and ensure all pins have locked in place (you will hear a clicking noise) once the pins have fully protruded through the holes.

Please ensure there is no sign of wear and damage.

Please ensure the user does not exceed the max user weight stated.

REF	Model No	Description	Handle Style	Handle to Armband	Handle to Ground	Max User Weight
091275965	121	Single Adjustable Crutches - Standard	Standard	230mm	660 - 940mm	160kg/25st
091356468	121T	Single Adjustable Crutches - Extra Long	Standard	230mm	830 - 1100mm	160kg/25st
091325323	121JDA	Double Adjustable Crutches - Childs	Standard	170 - 220mm	560 - 800mm	100kg/15¾st
091115799 (091534643)	121A (121A-D)	Double Adjustable Crutches - Standard	Standard	200 - 275mm	660 - 940mm	160kg/25st
091356443 (091525153)	121AT (121AT-D)	Double Adjustable Crutches - Extra Long	Standard	200 - 275mm	830 - 1100mm	160kg/25st
091356476	124	Single Adjustable Crutches - Standard	Ergonomic	230mm	690 - 1100mm	160kg/25st
091356484	124A	Double Adjustable Crutches - Standard	Ergonomic	200 - 275mm	690 - 940mm	160kg/25st
091356492	124AT	Double Adjustable Crutches - Extra Long	Ergonomic	200 - 275mm	830 - 1100mm	160kg/25st
091356385	119E	Fixed Height Crutches	Ergonomic	235 - 310mm	1030mm	190kg/30st
091439199	125	Single Adjustable Crutches - Standard	Comfy	230mm	660 - 940mm	160kg/25st
091531540	125T	Single Adjustable Crutches - Extra Long	Comfy	230mm	830 - 1100mm	160kg/25st
091439207	125A	Double Adjustable Crutches - Standard	Comfy	200 - 275mm	660 - 940mm	160kg/25st
091531557	125AT	Double Adjustable Crutches - Extra Long	Comfy	200 - 275mm	830 - 1100mm	160kg/25st

Length Adjustment

When the crutch length is correct the user should be able to maintain an upright posture with the elbow slightly flexed. In this way body weight is taken through the crutch by pushing down on the crutch when walking.

To adjust a crutch fitted with a pin clip, press the two buttons in. Push the leg in and turn slightly. Adjust the leg to the required position, then turn and relocate, ensuring buttons are protruding through the holes.

Where the affected leg can take some weight move the crutches one pace forward then place the affected leg on the ground with the crutches, which will support some of the body weight. Step through with the unaffected leg.

Getting up and sitting down

The crutch should be near the chair or held in one arm. Lean forward, push up on one crutch and the opposite arm of the chair to assist in standing. Once standing up ensure the arm is secure inside the cuff.

When returning to the chair, turn round, feel the chair behind the legs. Remove both arms from the crutches, hold them in one hand or place them near the chair. Place hand(s) gently on the chair arms, bend forward and lower gently into the chair.

Walking up and down stairs

Hold onto a handrail if at all possible. Hold the spare crutch together with the handle of the remaining crutch. The crutch goes on the same step as the affected leg. The unaffected leg should lead when going upstairs, and the affected leg leads when coming down.

Usage

Walking

Where the affected leg cannot take any weight keep it off the ground by holding it somewhat in front with the knee slightly bent with arms through cuffs and gripping handles. Move the crutches one pace forward, slightly apart and level then move the body forwards between the crutches allowing the unaffected leg to swing through.

Safety and Maintenance Information

Do not put any force on the cuff. The cuff is for stability when using the crutch and is not load bearing. It should not be used to lever or aid rising from a seated position or ascending/descending from a flight of stairs.

If the crutch is dropped or slips from the hand accidentally, please examine the plastic components (handle, cuffs) for any signs of cracks which may have occurred before using.

Please check the ferrules for wear and tear regularly. On a monthly basis please check screws on the handle are secure. If tightening is required, ensure the screws are not over-tightened.

All crutches are designed for minimum maintenance. Manufactured from anodised aluminium tubing they require the occasional wipe-over with a soft cloth. Do not use harsh abrasive materials!

If disinfection is required, clean as above and wipe over with a solution of 1,000 parts per million available chlorine disinfectant (haz tabs/precept), rinse and dry. Decontamination process that involves immersion is not recommended.


Do not store in sub-zero temperatures.

Weekly check for elongated or stressed adjustment holes, bent sections or split handles. Discard the product if any of these faults are found. Split, worn or loose ferrules must be replaced. Check for rust on the adjustment spring buttons.

Avoid wet floors, outdoor hazards such as wet leaves and ice. Wear flat supportive shoes. Do not exceed the maximum user weight stated on the product label.

Warranty Information

The Days range of adjustable crutches are warranted for one year from the date of purchase against faulty workmanship or materials. Please contact your supplier/dealer.

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These user instructions are available in large print, please contact Performance Health Customer Services for a copy. Please retain these instructions for future reference.

Symbols Glossary located at: <http://www.performancehealth.com/medical-device-glossary>

 Performance Health International Ltd.
Nunn Brook Road, Huthwaite,
Sutton-in-Ashfield, Nottinghamshire,
NG17 2HU, UK
UK:
Tel: 03448 730 035
www.performancehealth.co.uk
International:
Tel: +44 1623 448 706

 Performance Health France
13 rue André Pingat CS 10045
51 724 Reims Cedex, France

Distributed in the USA by:
Performance Health Supply, LLC
W68 N158 Evergreen Blvd.
Cedarburg, WI 53012 USA
+1 800-323-5547
www.performancehealth.com

Distributed in Australia by:
Ausmedic Australia, Trading as
Performance Health ANZ
Unit 3, 3 Basalt Road
Pemulwuy NSW, 2145, Australia
1300 473 422
www.performhealth.com.au

Distribué en France par:
Performance Health France
13 rue André Pingat CS 10045
51 724 Reims Cedex, France
00(33) 03 10 00 79 30
www.performancehealth.fr

EC REP
CEpartner4U
Esdoorlaan 13
3951 DB Maarn
The Netherlands
www.cepartner4u.com

**PERFORMANCE
HEALTH**