

**3/32" 38% Superperf**

**3/32" 19% Optiperf™**

**1/8" 1% Perf**

**1/16" 42% Perf**

**1/8" 42% Perf**

**1/8" 19% Optiperf™**

**1/16" 19% Optiperf™**

**1/16" 13% Ultraperf®**

## Choosing solid versus perforated material:

**Solid** for extra rigidity.

**Perforated** (1%, 2% or 21%) for slight ventilation and added comfort without compromising rigidity.

**Ultraperforated** (13%), **Optiperf™** (19%), and **Super perforated** (38% or 42%) for lighter weight, exceptional ventilation and greater comfort and compliance. For additional support, select a thicker material.

See below for details on our  
"Perforations Pack."



**Perforations Pack.** A great way to sample the different perforation patterns available from Rolyan. Each piece of Aquaplast®-T and Watercolors™ measures 6" x 9" (15cm x 23cm) and has a non-stick coating.

Thickness	Perforation	Color
1/8"	42% Superperf	White
1/8"	19% Optiperf	Spring Green
1/8"	1% Perf	White
3/32"	19% Optiperf	Purple
3/32"	38% Superperf	White
1/16"	19% Optiperf	Royal Blue
1/16"	13% Ultraperf	Lavender
1/16"	42% Superperf	White


## Why "buy a case"?

Our most popular sheet materials are available in case quantities to offer savings and convenience. We even offer mixed cases (pages 472, 474 and 483).




We offer a variety of Rolyan® thermoplastics with an antimicrobial component built in for odor reduction.


## Rolyan® Thermoplastics at a glance... Choosing the appropriate thickness:




- Half the weight of 1/8" splinting material
- Easy to cut when cold
- Ideal for pediatrics, finger splints and hand-based splints



- Lightweight with moderate support
- Minimizes bulkiness
- Perfect for a wide range of static and dynamic splints including arthritis splints and thumb, hand and forearm splints




- Traditional thickness for greater support
- Well-suited for upper and lower extremity splints, spasticity splints and fracture bracing



- Maximum rigidity for body jackets, spasticity splints, lower extremity splints and fracture braces

## Questions about splinting?

Call our Splinting Hotline at 1-800-228-3693, x 8771 or email us at [splinting@patterson-medical.com](mailto:splinting@patterson-medical.com)



Have an idea for a new product?  
Need something special? Just let us know  
and we'll be happy to work with you!