









Our most popular walking frames are now available in vibrant metallic Red

These new RED Days Walking Frames were developed to meet demand from the NHS market for use with patients with dementia.

Bright colours can highlight important objects, making them easier for people with dementia to find and use making this range ideal to help increase the mobility for patients with dementia and reduce the risk of potential falls.

	Wheeled and Non-Wheeled versions available		Hard anodised finish for durability
	Comfortable latex free PVC handgrips		Low maintenance and easy to wipe clean
	Lightweight aluminium frame		Simple push button height adjustment
	Non-marking rubber ferrules		Maximum user weight <b>160 kg</b> / <b>25 st</b>


**Days Adjustable Height Walking Frames** 

**1 Domestic/Narrow** 


	<b>Size</b>	
09 155 3890	Small	£54.62   3   No
09 155 3866	Medium	£418.88   19   No
09 155 3882	Large	£56.49   3   No

**2 Standard Hospital** 


	<b>Size</b>	
09 155 3775	Small	£895.47   20   No
09 155 3767	Medium	£380.6   14   No
09 155 3759	Large	£91.8   6   No

**3 Domestic/Narrow Wheeled** 

	<b>Size</b>	
09 155 3874	Small	£326.34   18   No
09 155 3841	Medium	£3066.19   111   No
09 155 3858	Large	£598.72   30   No

**4 Standard Hospital Wheeled** 

	<b>Size</b>	
09 155 3809	Small	£663   24   No
09 155 3791	Medium	£3209.26   132   No
09 155 3783	Large	£1143.53   42   No

**5 Spares For Above (Not Illustrated)** 

09 135 6120	25mm PVC ferrule	£3420.68   6707   No
09 155 3981	Wheeled Extension Legs	£0   0   No



Specifications						
	Size	Overall Width B	Width at Front	Depth C	Min. Height A	Max Height A
09 155 3890	Small	59cm (23½")	41cm (16")	46cm (18")	70 cm (27½")	78 cm (30¾")
09 155 3866	Medium	60cm (23¾")	42cm (16½")	46cm (18")	80.5 cm (31¾")	88.5 cm (34¾")
09 155 3882	Large	62cm (24½")	42cm (16½")	46cm (18")	87.5 cm (34½")	95.5 cm (37½")
09 155 3775	Small	65cm (25½")	46cm (18")	49cm (19¼")	67cm (26½")	75cm (29½")
09 155 3767	Medium	66cm (26")	47cm (18½")	55cm (21½")	77cm (30¼")	85cm (32½")
09 155 3759	Large	67cm (26½")	51cm (20")	57cm (22½")	87cm (34¼")	95cm (37½")
09 155 3874	Small	59cm (23½")	50cm (19¾")	49cm (19¼")	73cm (28¾")	81cm (31¾")
09 155 3841	Medium	60cm (23¾")	51cm (20")	49cm (19¼")	83.5cm (33")	91.5cm (36")