

Vinco®-MedWipe

Healthcare Biodegradable & Plant Based Plastic Free Wet Wipes For Cleaning & Disinfection



- ✓ Effective and safer on Hands, Surfaces & Equipment
- ✓ Provides Broad Spectrum Cleaning & Disinfection
- ✓ Unique Green Biodegradable and Plant-based Formula
- ✓ Kills Viruses & Bacteria Fast



Biodegradable



NEW AND
IMPROVED WIPE
RANGE FOR
HEALTHCARE



CALL NOW
03448 730 035



VISIT US AT
PERFORMANCEHEALTH.CO.UK

Vinco®-MedWipe Healthcare Biodegradable & Plant Based Plastic Free Wet Wipes For Cleaning & Disinfection

Anti-viral and anti-bacterial wet wipes for use on hard surfaces and non-invasive medical devices.

Introducing you to the Vinco-MedWipe with a unique formulation. Most conventional products have formulations that are based on actives such as QUAT, PHMB, Amine, Triamine, and Chlorine which can all aid antimicrobial resistance. Testing has also shown that these only provide limited efficacy and spectrum for disinfection where killed or injured organisms are still viable to potentially mutate to superbugs or transfer their resistance elsewhere.

Some of these actives can also be harmful to humans, animals and the environment as well as being corrosive to skin, surfaces and equipment. Using the Vinco®-MedWipe with a unique disinfection formula you can be sure that you can carry out safe cleaning and disinfection, also knowing the MedWipe safeguards users from today's infection risks.

Vinco®-MedWipe also prevents mutation and reduces antimicrobial resistance [AMR] and protects our antibiotics as well as human, animal & environmental health.

CONVENTIONAL PRODUCTS

Soap

Low - level disinfectants PHMB, Quats, Triamine etc.

Ethanol

Chlorine

Vinco®-MedWipe

CLEAN & DISINFECT WITH Vinco-MedWipe to REPLACE 97% of conventional products



Key benefits of using Vinco®-MedWipe:

- // Compliant with over 20 British standards.
- // Approved by GB HSE under the Simplified Authorisation.
- // Safeguards users and personnel from today's infection risks.
- // Dermatologically Tested
- // Prevents mutation and reduces antimicrobial resistance [AMR].
- // Biodegradable plastic free range
- // Non-hazardous solution, no rinse and no PPE required.
- // Kind to young and sensitive skin, vegans and the environment.
- // Made In Britain
- // Plant-based product - no hazard warnings required.
- // Protects our antibiotics as well as human, animal & environmental health.

For more information on the Vinco®-MedWipe range: Visit: www.performancehealth.co.uk
| Call: 03448 730 035 | Email: ukcustomersupport@performancehealth.com



7028217
 Vinco-Medwipe Tub
 Cleaning & Disinfection
 Bio Wipes
 200 sheet | 20x20cm



7028218
 Vinco-Medwipe Bucket
 Cleaning & Disinfection
 Bio Wipes
 500 sheet | 20x20cm



7028220
 Vinco-Medwipe Refill Bag
 Cleaning & Disinfection
 Bio Wipes
 500 sheet | 20x20cm



7028219
 Vinco-Medwipe Pouch
 Cleaning & Disinfection
 Bio Wipes
 200 sheet | 20x20cm



7028215
 Vinco-Medwipe
 Tub Cleaning &
 Disinfection PP Wipes
 200 sheet | 20x20cm



7028216
 Vinco-Medwipe
 Pouch Cleaning &
 Disinfection PP Wipes
 200 sheet | 20x20cm

- ✓ Safe on hands and surfaces
- ✓ Plant-based formulation
- ✓ Plastic free biodegradable wipes
- ✓ Kills 99.999% of all bacteria in 1 minute
- ✓ Kills 99.99% of all viruses in 2 minutes including Coronavirus
- ✓ Biodegradable plastic free range
- ✓ Kills yeast in 5 minutes
- ✓ Tested to BSEN standards
- ✓ Free from PHMB & QUAT
- ✓ Rinse free

Code	Description	Packaging Format	Sheet Size	Qty. Per Pallet
Vinco-MedWipe Bio Range				
7028217	Vinco-MedWipe Cleaning & Disinfection Bio Wipes 200sheet	TUB	20x20	432
7028218	Vinco-MedWipe Cleaning & Disinfection Bio Wipes 500sheet	TUB	20x20	192
7028220	Vinco-MedWipe Cleaning & Disinfection Bio Wipes 500sheet	EQP	20x20	192
7028219	Vinco-MedWipe Cleaning & Disinfection Bio Wipes 200sheet	EQP	20x20	432
Vinco-MedWipe PP Range				
7028215	Vinco-MedWipe Cleaning & Disinfection PP Wipes 200sheet White	TUB	20x20	423
7028216	Vinco-MedWipe Cleaning & Disinfection PP Wipes 200sheet White	EQP	20x20	432

Suitable for:

The Vinco®-MedWipe range is suitable for use on hard surfaces and non-invasive medical devices. Designed and developed from science it emphasises cleaning and disinfection. This multi-purpose wipes range is formulated from a plant-based active and coupled with a high performing 100% biodegradable viscose wipe. The plant-based formulation is free from harmful biocides such as QUAT, PHMB, Amine, Triamine, and Chlorine which are all widely used across other conventional products. This wipe is powerful, safe and delivers broadspectrum disinfection of multi-resistant bacteria, viruses, yeast and mycobacteria approved for healthcare and the dental sectors.

Efficacy	Product Type Group	Contact Time
Bactericidal Efficacy		
EN1276	2 & 4	1 min
EN13697	2 & 4	1 min
EN13727	2	1 min
EN14561	2	1 min
EN16615	2 & 4	1 min
EN1500	1	1 min
PAS2424 - Bacteria 24-hour residual	2 & 4	1 min
Virucidal Efficacy		
EN16777 - enveloped and non-enveloped	2 & 4	2 min
EN16777 - 24hr residual enveloped	2 & 4	2 min
EN14476	2 & 4	2 min
Yeasticidal Efficacy		
EN1650	2 & 4	5 min
EN13624	2	3 min
EN14562	2	5 min
EN13697	2 & 4	5 min
EN16615	2 & 4	1 min
PAS2424 - Yeast 24-hour residual	2 & 4	5 min
Mycobactericidal Efficacy		
EN14348	2 & 4	10 min
EN14563	2	10 min